Use of Environmental DNA to Detect and Quantify Brook Trout Populations in Adirondack Mountain Streams

Barry P. Baldigo, U.S. Geological Survey Lee Ann Sporn, Paul Smith's College Scott D. George, U.S. Geological Survey








Background

- Fish communities were surveyed 40 SW Adirondack streams by the USGS and NYSDEC, summer 2014-15
- Density and biomass of Brook Trout populations were quantified
- eDNA samples were collected in late summer 2015








Base from National Geographic, Esri, DeLorme, HERE, UNEP-WCMC, USGS, NASA,

Primary Objective

- Determine if eDNA can be an effective tool for assessing the distribution and abundance of Brook Trout in remote/headwater streams of the region.
- Related research goals were to:
 - a. Compare the effectiveness of water and sediment to retain Brook Trout DNA
 - b. Evaluate eDNA water-filtration volumes
 - c. Assess the accuracy of eDNA to predict presence/absence
 - d. Assess the accuracy of eDNA to predict density and biomass
 - e. Weigh costs and benefits of standard e-fish survey vs. eDNA sampling efforts



Fish surveys

- Selected study streams that represented a wide range in acid-base chemistry
- Isolated sampling reach using blocking seines
- Collected all fish using 3+ electroshocking passes
- Identified & measured lengths and weights of each fish
- Collected basic habitat data at each study reach
- Computed metrics for fish communities and Brook Trout populations






Brook Trout eDNA

- Decontaminate all gear with 10% bleach & rinse well with site water
- Collect 2-300 g sediment, filter up to 6L water thru 1.5 μ m GFF, and freeze
- Extract, purify eDNA from filters and sediment using DNA isolation kits




- Detect copies of Brook Trout DNA with the TaqMan® MGB assay using a Step OneTM qPCR system, Environmental Master Mix 2.0, and reaction components: 250 nM BRK2 probe; 900 nM BRK2 primers; VIC -labeled internal positive control (IPC) primer/probe; and IPC target DNA with forward primer = ccacagtgcttcaccttctatttcta, reverse primer = gccaagtaatatagctacaaaacctaatagatc, and FAM-labeled probe = actccgacgctgacaa
- Each assay used 1μ I of eluted DNA and all samples were run in triplicate
- Average cycle threshold (Ct) value for each sample was calculated from all triplicates with < 40 amplification cycles
- The relative concentration of Brook Trout DNA (ng/ μ l) in each sample was estimated from a standard curve ($10^{(Ct-26.685)/-3.615}$; $R^2 = 0.998$) between mean Ct and total genomic DNA (purified from local Brook Trout tissues)

RESULTS



RESULTS



eDNA RESULTS

Predicting presence:


- After the first survey, eDNA correctly detected Brook Trout at 25 of the 30 sites where populations were observed (83% classified correctly)
- After resampling 2 of 5 misclassified sites, eDNA detected Brook Trout at 27 out of 30 sites (90% classified correctly)
- The remaining three sites were not resampled

Predicting absence:

- After the first survey, eDNA did not detect Brook Trout at 9 of 10 sites (90% classified correctly) where populations were not found
- After resampling 1 misclassified site, eDNA did not detect Brook Trout at all 10 sites (100% classified correctly)
- Overall, eDNA correctly classified the P/A of Brook Trout populations in 93% of study sites

RESULTS

Ability to predict Brook Trout density and biomass



Discussion - Conclusions

- Sediments did not retain Brook Trout DNA in headwater Adirondack streams
- Decontamination & rinsing procedures for eDNA samples need to be thorough
- Large sample volumes and/or field duplicates are crucial for detecting rare targets (threshold)
- Costs to analyze an eDNA sample (\$10-50) can be a fraction of that for a crew to survey fisheries (\$0.5-3K)
- eDNA is an effective tool for detecting the P/A and relative abundance of Brook Trout populations in small headwater streams
- eDNA is appropriate for large/broad/synoptic surveys; standard surveys are generally more effective when characterizing the health of entire fish assemblages



QUESTIONS?

bbaldigo@usgs.gov

Special thanks:

Jason Ball, Noel Deyette, Luis Rodriguez, Nick McCloskey, Howard Simonin, Eric Paul, Ryan Collins,

and Spencer Bruce





Department of Environmental Conservation



Paul Smith's College